+

Lumi

EN	Light therapy mask	2
FI	Valohoitonaamio kasvoille	9
SE	Ljusterapimask	16
LT	Šviesos terapijos kaukė	23
LV	Gaismas terapijas maska	30
EE	Valgusteraapia mask	37

EN Important information

Users of this device can obtain benefit from using it only if they read and comprehend all information in the user manual. You must follow all the instructions in the user manual. Please read the accompanying user manual thoroughly and follow the directions stated in the manual to assure maximum safety during operation.

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In the package you will find

How Light therapy mask works

Who should not use the Light therapy mask? Contraindications

General conditions

Warning

Precautions

How to use

Before starting to use

Skin test

Recommended treatment schedule

Using the device

Colour modes

After use

Common skin reactions

Cleaning

Aftercare

Recycling

Technical specifications

Troubleshooting

Warranty and support

Welcome

Welcome to the beauty world of Lumi! Lumi uses non-invasive LED light therapy, known as one of the most effective methods to help promote skin rejuvenation, reduce the appearance of pigmentation, wrinkles and acne. The Light therapy mask is designed to deliver light directly to the skin, ensuring maximum exposure with a choice of light in 3 different colors plus the infrared region of the light spectrum for the targeted action you want.

In the package you will find

- Light therapy mask;
- Control panel;
- Power cord;
- Silicone eye protection pads;
- 3x adjustable straps.

How Light therapy mask works

Photonic light waves of different lengths reach the deepest layers of the skin, blood vessels, and lymphatic system. Depending on the chosen program, light energy is converted into cellular energy in a specific layer of the skin. This conversion initiates a number of metabolic actions, such as the formation of new capillaries, the stimulation of collagen production, and the production of adenosine triphosphate (ATP). The main function of ATP is to supply our cells with useful energy.

Who should not use the Light therapy mask? Contraindications

General conditions

- Never use the Light therapy mask if you are pregnant, planning on becoming pregnant or breastfeeding. The Light therapy mask has not been tested on pregnant or breast-feeding women and therefore the risk to the fetus, neonate or pregnant women is unknown.
- Never use the Light therapy mask to treat any other conditions apart from those listed in the indications for use. The Light therapy mask has not been tested for any other conditions than those listed, and the risk is unknown.
- Never use the Light therapy mask if you suffer from Lupus erythematosus, photosensitive eczema or Albinism. If you use the Light therapy mask to treat Lupus erythematosus, photosensitive eczema or Albinism you may cause a severe skin reaction.
- Never use the Light therapy mask if you suffer from any photosensitive disorder (sensitization to light). If you use the Light therapy mask and you suffer from a photosensitive disorder you may cause a severe skin reaction.
- Never use the Light therapy mask if you are taking any medication that can cause photosensitivity. If you use the Light therapy mask and you are taking any medication that can cause photosensitivity you may cause a severe skin reaction.

Note: Photosensitivity is a common side effect of various medications. These can include certain antibiotics, chemotherapy drugs, and diuretics. If you are unsure about any medication you may be taking, consult your healthcare provider.

Note: We recommend that you seek medical advice before using the device if you have any doubts or concerns.

Warning

- DO NOT use the Light therapy mask in any way that is not exactly as specified in these instructions for use. The use of the Light therapy mask other than specified in these instructions for use has not been tested and the risk is unknown.
- DO NOT use a Light therapy mask without first removing creams, serums, lotions, and make-up BEFORE the procedure. Some ingredients in these products may cause light sensitization or may reduce the effect of the treatment. If you use the Light therapy mask without first removing these products, you may cause a skin reaction, damage your device, or you may reduce the effect of the treatment.
- TURN OFF and REMOVE the mask immediately if it feels hot or your skin feels uncomfortably warm.
- NEVER use the near water or where water can be split into the equipment. DO NOT use the mask if you think it has come into contact with water.
- Keep the Light therapy mask away from HOT surfaces.
- DO NOT attempt to modify the Light therapy mask in any way. Modifying the Light therapy mask will void your warranty. Modifying the mask could cause an electric shock when the system is used again. The system contains no serviceable parts.
- DO NOT use the Light therapy mask if the mask, controller or cord is damaged in any way.

Note: Stop using the Light therapy mask if you encounter any prolonged skin redness or discomfort that lasts for more than 2 hours. If a problem occurs for more than 2 hours you might be light sensitive or have another condition. We recommend you consult with a medical professional for advice.

Precautions

- Before each use, make sure the device is switched off and clean the Light therapy mask. All surfaces of the mask may be cleaned by following the cleaning instructions.
- Do not use any other cleansing products other than those stated in the cleaning instructions. Use of other cleansing products may damage the mask.
- Keep the Light therapy mask at room temperature in dry conditions, away from direct sunlight.

How to use

Before starting to use

1. Fully charge the device:

- Connect the controller to the power source using the provided charging cable.
- Charge the controller separately with the wall socket cable. Charging while connected to the device may damage the battery.
- The 3 LEDs in the control unit will start to flash to indicate charging and remain steady when fully charged. Charging typically takes about 4 hours.
- Ensure the device is fully charged before the first use for optimal performance.

2. Check all components:

- Refer to the section "In the package you will find" to ensure all listed parts are included and undamaged.
- If any parts are missing or damaged, do not use the device and contact your provider for assistance.

Skin test

You must follow the instructions below before your first treatment session to see if you are light sensitive:

- Secure the Light therapy mask to your face, allowing the LEDs to deliver light directly to your skin.
- Turn the Light therapy mask on (see the "Handling the device" section of the instructions) – and wear it for 5 minutes. At any time if your Light therapy mask feels HOT or uncomfortable STOP using the device immediately.
- Wait for 24 hours.
- If you experience pain, discomfort or redness that lasts for more than 2 hours you might be light sensitive. Please consult with a medical professional for advice.

Recommended treatment schedule

- For the best results we recommend 10-25 minute treatments, 3-4 times per week (with 2 days rest) for a minimum of 12 weeks.
- Maintenance treatments (1-2 treatments per week) can then be carried out as required once the desired results are achieved.
- **Reminder:** For enhanced results, the **Lumi face mask and Lumi neck and décolleté mask can be used simultaneously**. Connect both device cables to the controller's dual ports to experience the combined benefits of dual treatments in the comfort of your home.

Using the device

STEP 1

Prepare your face and connect the mask.

Gently remove makeup or cleanse your face and allow it to dry naturally. Once your skin is ready, connect the Light Therapy Mask to the control panel.

STEP 2

Attach the silicone eye pads.

Place the silicone eye pads into the mask's eye openings, ensuring the narrower ends point outward (as shown in the image).

STEP 3

Secure the mask using adjustable straps.

Thread the straps through the slots on the mask, adjust their length, and fasten them around your head for a snug but comfortable fit.

STEP 4

Turn on the mask and choose settings.

Press and hold the power button (~2.5 seconds) to turn on the device.

- Select a light color using the M button: Red (R), Blue (B), Yellow (Y) or Combined (G).
- Adjust light intensity with the E button, starting with a low level and increasing gradually.
- Set the session duration (recommended 10–25 minutes) using the T+ and T- buttons.
- Activate near-infrared light (NIR) by pressing the NIR button if needed.

STEP 5

Wear the mask and begin the session.

Put on the mask and ensure it is securely in place. Relax during the procedure and avoid looking directly at the light.

STEP 6

Clean the mask and follow up with skincare.

After the session, remove the device and wipe the device clean (see cleaning instructions). Apply your regular skincare products, including sunscreen if the session was in the morning.

STEP 7

Repeat as needed.

Use the mask 3–4 times a week for 10–25 minutes. Once results are visible, reduce the frequency to 1–2 sessions per week.

Colour modes

Three photonic light programs allow you to select the desired effect depending on the condition of your skin:

- *Red (smoothes wrinkles and scars, tightens sagging skin).* 630 nm

The widest spectrum of waves stimulates blood circulation and collagen production. Helps retain moisture in the skin, tightens and makes the skin more elastic.

- *Blue (fights acne and facial rashes. Antibacterial, anti-inflammatory).* 460 nm

The program with antibacterial properties for inflammatory acne is especially helpful in fighting the acne-causing bacterium *P. acnes*. Stabilises sebum production.

- *Yellow (reduces redness of the skin, dark circles under the eyes).* 590 nm

Detoxifies the skin, stimulates the lymphatic system, effectively constricts dilated capillaries. Fights red, rosacea, erythema.

- *Near infrared light (improves the general condition of the skin).* 850 nm

Stimulates cell renewal, improves blood circulation in capillaries. Increases the level of moisture in the skin.

After use

Common skin reactions

Immediately after treatment, skin may appear red. Skin should return to normal within 1-2 hours.

Cleaning

Make sure the Light therapy mask is unplugged and turned off. Then use a damp, soft cloth with a mild alcohol-free cleaning detergent to clean the facial area. Repeat until the mask is clean and then dry with a clean dry cloth. Do not use alcohol or harsh chemicals or any acidic cleaning materials to clean the Light therapy mask as these may damage the mask and affect how it works.

Aftercare

- Store away from high temperatures and moisture. Avoid storing the Light therapy mask in dusty conditions, in direct sunlight or in damp and humid conditions.
- Store the Light therapy mask at room temperature in dry conditions.
- Always adjust your Light therapy mask using the head straps, never the cord that connects the mask to the controller.
- Handle with care to avoid accidental damage.
- Avoid tugging, pulling or bending the cord that connects the Light therapy mask to the controller.

Note: After cleaning, we recommend you to store the product in its original box.

Recycling

Follow your country's rules for the separate collection of electrical and electronic products and rechargeable batteries. Correct disposal helps prevent negative consequences for the environment and human health.

Technical specifications

- Device size: 327 * 210 * 5.8 mm;
- Package size: 341 * 230 * 49 mm;
- Weight of the device: 1,2 kg;
- Voltage: DC 5V
- Power Supply: USB-powered (via external USB adapter or computer USB port)
- Amount of light emitted: 3;
- Number of LEDs: 280;
- Storage conditions temperature: -5°C to 55°C;
- Storage conditions humidity: 30% to 95% without condensation.

Troubleshooting

Problem	Solution
The Light therapy mask does not turn ON	Please check if the controller is charged; that the mask is connected to the controller properly.
The Light therapy mask is flicking	Please check if the controller is charged; that the mask is connected to the controller properly.
The Light therapy mask mask is hot	This is normal, the face mask becomes hot during a treatment that is designed to open the pores of the skin. If the face mask becomes too hot, please turn off and wait for it to cool down.

Note: If the actions suggested fail to fix the problem, or the problem reoccurs, and your device is within the warranty period, please return to your provider. Do not attempt to investigate any further or modify the device in any way. The Light therapy mask has NO parts you can repair yourself.

This product is not intended to diagnose, cure or prevent disease. The Light therapy mask makes no claims, representations or otherwise warranties the ability of this product to cure any physical, skin or mental conditions. We do not claim that the Light therapy mask prevents disease or cures ailments or conditions. Our products must be treated as fitness products or general wellness products.

Manufacturer:

Shimeng Industrial (Shenzhen) Co., LTD
Floor 2, Building 9, ZhengZhong Science Park,
QiaoTou Community, Fuyong Street, Baoan, Shenzhen
Tel: 0755-26514485
Fax: 0755-26514485-802
Email: sales3@freeflying.cn

Warranty and support

This sheet contains information about the warranty conditions of the product. For detailed information and instructions on the use and maintenance of the device, see the instructions included with the order.

Warranty claim form

Product name	
Date of purchase	
Business name, address	

1. Warranty period

The warranty is valid from the date the purchaser purchased this product from the seller. The warranty period is 24 months.

2. Warranty conditions

The warranty is valid upon presentation of the original product purchase document.

If the product is defective during the warranty period for any reason not covered by the warranty, the manufacturer undertakes to repair or replace the product. The warranty does not cover regular inspection, maintenance and repair or replacement of normally worn parts.

The warranty does not cover damage or malfunctions caused by:

- Improper use of the product and / or procedures not in accordance with the manufacturer's instructions for use, technical requirements or safety standards, including improper maintenance, falling to the ground and severe shocks.
- corrosion, dirt, water.
- Repairs performed outside of warranty service.
- Improper cleaning of the appliance, notwithstanding the recommendations given.
- Accidents or incidents or any causes beyond the control of the manufacturer.
- If the product does not meet the standards or specifications of a specific country other than the country of purchase.

Any attempt to adapt the product to the technical or safety requirements of other countries will void all warranty service. If no defects are found in the product, all service-related costs are paid by the buyer.

The warranty claim must be submitted by email within two weeks of the defect occurring, describing in detail the basis for the warranty claim. In the event that the application is recognized as valid, the guarantor may, at its discretion:

- replace the defective product with a new or refurbished product,
- free of charge to find any defects found,
- to pay financial compensation when the elimination of defects will not be possible or if the elimination of the defect requires excessive costs depending on the value of the defect-free product and the type or significance of the defect found.

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